



MEDSTONE

Zogear[®] Child's Face Mask

Child mask is fluid resistant and provides the child protection against large droplets, splashes or sprays of bodily or other hazardous fluids. Protects friends, family and others from the child's respiratory emissions.

Specifications:

- PFE: > 99%
- BFE: > 99%
- Breathability (Delta P): < 5.0 H2O/cm²
- Flammability: Class I

Certifications:

- CE
- FDA
- ISO13485

Characteristics:

- 100% Latex-free soft stretchable ear-loops
- Fluid resistant inner layer
- Adjustable nose piece
- High-standard breathability, fiberglass free, hypoallergenic
- Comfortable inner lining, designed for sensitive skin



MEDSTONE

Medstone
Amstelvein 28A
Mondriaan Tower (13th Floor)
1096 BC Amsterdam
The Netherlands

Company website: www.medstone.nl
Online shop: shop.medstone.nl
Contact: info@medstone.nl

Disclaimer: No rights can be derived from the information provided on this sheet
© 2020 All Rights Reserved.
Medstone[®] is a registered trade name of Fresco Farma B.V.